

GKF 2000/2500 Cleaning Station Upgrade Kit



BOSCH



P/N 8111338599

STICKY, DUSTY PRODUCTS CUTTING DOWN YOUR YIELDS?

Bosch offers a new Cleaning Station Upgrade Kit to intensify the segment cleaning ability of your GKF 2000/2500. Your current plumbing and timing are modified to maximize the cleaning ability of your GKF 2000/2500.

- Longer cleaning pins allow both lower and upper segments to be thoroughly cleaned
- Minimized capsule separation problems due to less product residue left on segments
- Easy on-site installation
- Increased yields by increasing capsule separation efficiency

Bosch Packaging Technology
8700 Wyoming Ave. N.
Minneapolis, MN 55445
763-424-4700

www.boschpackaging.com

