

TL Large 8 Pump Filler Preventative Maintenance Kits



BOSCH

At Bosch Packaging Services, we understand the importance of a good preventative maintenance program.


Such a program requires the routine replacement of certain wear parts. We have identified select bearing, seals and wear parts which we recommend be replaced in a preventative maintenance program for the TL large 8 pump fillers. Utilization of our kits will help ensure the long term fill accuracy and reliability of your filler.

Note: These kits do not include the cams since they are different for every machine. When requesting a quote, please request cams as these are a critical part of your filler.

Service is recommended but not required.


All of the minor preventative kit parts are included in the major kit.

Each minor or major kit includes enough parts to rebuild (1) pump mounting assembly. Additional pump mounting assembly rebuild kits are available at additional cost. Please call for part numbers and pricing.



* 76600 Minor preventative maintenance kit recommended yearly or every 3,000 hours of use (estimated one day or 8 hours of service).

* 76601 Major preventative maintenance kit recommended every 3 years or 6,000 hours of use (estimated five days or 40 hours of service).



Major Rebuild kit includes all items pictured. Minor Rebuild kit only includes those pictured on the right.

